

RECONSTRUCTION *refined*



About AeroForm

AeroForm is designed for women who want to take their breast cancer head on, and want to undergo breast reconstruction on their terms and their timeline. We're committed to helping women regain the control that has been taken away during their journey.

AeroForm is a revolutionary **needle-free**, breast tissue expander for women who elect to have reconstructive surgery following a mastectomy. Unlike traditional saline expanders, AeroForm is designed to bring control, convenience and comfort to the breast reconstruction process.

CONTROL

AeroForm is the only tissue expander that empowers **you** to control your expansion.

CONVENIENCE

With AeroForm, **you** can expand whenever and wherever you want.

COMFORT

AeroForm gives **you** the first and only needle free expander for breast reconstruction.

How AeroForm works

During surgery, your surgeon will place the AeroForm Expander under the skin in the space created by the mastectomy procedure. You will be able to start the expansion process once you have healed under the direction of your doctor. With AeroForm, you will be given your own personal Dosage Controller so that you can comfortably dose yourself up to three times a day. And of course, AeroForm is designed to give you the natural shape you are looking for.




“I felt empowered and in control; cancer can make you feel defeated and little things like self-dosing helped me feel more in control.”

—AeroForm Clinical Trial Patient


AeroForm data

AeroForm was tested in five separate clinical trials with over 250 women. The trials demonstrated:



PACE 1: n=7, PACE 2: n=33, ASPIRE: n=21, XPAND: n=130, XPAND II: n=50

The AeroForm difference



We are proud to bring you the latest technology in breast reconstruction. Ask your doctor if AeroForm is right for you.

* Ascherman, JA, et al. CO2-based vs. Saline Tissue Expansion for Breast Reconstruction: Results of the XPAND Prospective, Randomized Clinical Trial. Plastic and Reconstructive Surgery. December 2016 - Volume 138, Number 6. Consult with your physician if the AeroForm Tissue Expander System is right for you. Refer to product labeling information for potential adverse events, contraindications, warnings and precautions. The AeroForm Tissue Expander system is used for soft tissue expansion in breast reconstruction following mastectomy, for the treatment of underdeveloped breasts, and for the treatment of soft tissue deformities in the breast. The AeroForm Expander is intended for temporary subcutaneous or submuscular implantation and is not intended for use beyond six months.

© 2017. All rights reserved. AeroForm Tissue Expander System, the AeroForm logo, AirXpanders and the AirXpanders logo are registered trademarks and MMV is a trademark of AirXpanders, Inc. The products described in this document may be covered by one or more of the following U.S. Patents: 9,526,584; 8,394,118; 8,617,198 and 8,808,322. COL-0076 Rev A